



GLAMOURESSENCE™
"EMPOWERING BEAUTY. ENHANCING GLAMOUR"

2 IN 1

PICOSECOND + LONG PULSE

Tattoo Removal and Hair Removal Machine

532NM 755NM 1064NM 1320NM



Before and After Treatment Comparison

ENERGY:

3600MJ



3000W
Output Power



450PS
Pulse width

4

FOUR
Wavelengths

**NEW
UPGRADE**
14-inch screen



ADVANTAGES



LASER ARM DISPLAY

South Korea Imported 7 Joints
Light Guiding Arm



**TEC COOLING ENSURE
MACHINE CAN WORK
LONGER TIME**

IMPORTED LASER ROD

Two Standards (8) (7+8)
available




Heraeus



**HERAEUS LAMPS
GERMAN IMPORTED LIGHT SOURCE**



**ENERGY
3600MJ**



4 MODELS

- 1. Picosecond
- 2. Long Pulse
- 3. PTP
- 4. Multi Pulse



**NEW INTERFACE
ENERGY UPGRADE
CONFIGURATION
UPGRADE**

MODE SELECTION

4

MODES AVAILABLE



Picosecond

Tattoo removal treatment



Long Pulse

Hair removal treatment



PTP


Safer and more effective treatment of chloasma



Multi Pulse

Carbon peel treatment for skin rejuvenation and deep cleaning

LASER ARM DISPLAY



360° Rotation, Easy to Operate.

High energy output with low energy loss.

Maximum beam transmittance (**>95**).

7-section light guide arm imported from Korea.

Unlimited shots, full life span.

PROCESSOR HEAD




HONEYCOMB 755NM

Break The Melanin, More Dispersed! A Revolution Whitening.



1320NM

For Skin Rejuvenation, Shrink Pores, Skin Whitening.



RECTANGLES 755NM

Break The Melanin, More Dispersed! A Revolution Whitening.



1064NM / 532NM

For Dark Color Tattoo And Pigmentation Removal.



Powerful pigmentation and tattoo removal program.

ROD AND LAMP

Real




Φ 8 mm Laser Rod Made



Heraeus Lamps

German imported light source



Stable and reliable light output

Imported German laser lamp

Compact structure, Small laser volume

CAPACITOR UPGRADE




4 CAPACITORS


450V22000UF

SCREEN



14 INCH

-  Convenient operation
-  Fast response time
-  High-resolution screen display

BEFORE AND AFTER




- ✓ Not Only Pigment Removal
- ✓ But Also Hair Removal





COLOR SELECTION



Two colors available

 **black**  **white**

CUSTOMER REVIEW




Machine Parameters



NAME Picosecond Laser Tattoo Removal System	OEM&ODM Available
TYPE OF LASER Picosecond	OUTPUT POWER 3000w
WAVELENGTHS 532nm/1064nm,755nm honeycomb probe, 755nm Rectangular probe, 1320nm.	SCREEN 14 Inch Color Touch Screen
ENERGY 100-3600MJ (Adjustble)	SPOT AREA 2mm-10mm (Adjustble)
FREQUENCY 1HZ-10HZ (Adjustble)	VOLTAGE 110V/220V
COOLING SYSTEM Water Cooling + Air Cooling + TEC Cooling	PACKING Aluminum Alloy Box
WARRANTY 2 Years Warranty Free Technical Service	


WORKING PRINCIPLE

Picosecond laser




The PicoSecond Laser uses ultra-short pulses to hit the melanin with a great pressure, the melanin shatters into tiny dust-like particles. Because the particles are so small, they are more readily absorbed and eliminated by the body.

Nanosecond Technology VS Picosecond Technology



Traditional Laser

Traditional nanosecond lasers predominantly rely on photothermal action, delivering heat to the pigment and surrounding tissue.



Picosecond Laser

Picosecond lasers take advantage of pressure wave technology to shatter the target ink into tiny particles that are easily eliminated by the body.