

● 755+808+1064NM

DIODE LASER

With New AI Skin Analyzer

Hair Removal Machine



1600W
15*40MM



1200W
15*26MM






ADVANTAGES



1600W
Power



New AI 3D
Skin Analysis Technology



4K 15.6 Inch
Screen




Coherent
Bar



Tectonic
Design



Three
Wavelength




USB
Security Locks & Upgrades



Android
System



100
Million Shots



AI Skin Analyzer

Intelligent generate treatment parameter
Automatic skin & hair condition analysis

2 Million Pixels

The device supports the function of skin analysis before and after procedures, which support you provide more professional hair removal service.





GLAMOURESSENCE™
"EMPOWERING BEAUTY. ENHANCING GLAMOUR"



1600W Power 15*40MM



High power **1800W+1600W**



Handle Shots **100 Million Shots**





1200W Power

15*26MM

BEST COOLING SYSTEM

can reach to **-42°C**, won't burn skin.



Facial Tips



Diameter 8mm

Easily reach the parts that are hard to reach, including the ears, nostrils and eyebrows.

Target area:



Ear



Nose



Mouth Lips



Bikini



Upgrades & New Features



**Support storage
and query of treatment records**




**Support remote
control for renting**

The customer management system allows you to store customer records and exchange data between different clients.



GLAMOURESSENCE™
"EMPOWERING BEAUTY, ENHANCING GLAMOUR"

COHERENT BAR



16 Bar



12 Bar



COHERENT®


Our USA Coherent Bar

● More laser shots **100 millions**

● Real 3 wavelength **755+808+1064**

● High power each bar **100W**

● Stable quality longer lifetime




USB Internal Safety Lock & Upgrade



One-click push treatment parameters

Start your pocket beauty business from here

MOBILE INTERNET SYSTEM



Remote control system
Anywhere, Anytime, Remote control




User data management system
Simple process with clear parameters

4K 15.6 INCH SCREEN



16 Languages Options

Customizable welcome page and screen logo



Skin testing



User data

Android System

Latest wireless protocols

Latest Bluetooth 5.0


HANDLE & SCREEN OPERATION

Synchronization of handle operation data with the screen





BEFORE & AFTER






MACHINE PARAMETERS

LASERWAVELENGTH	755nm+808nm+1064nm
SCREEN	15.6 inch
HANDLE LASER POWER	1200W+1600W
PULSE DURATION	1-400ms
ENERGY	1-150J/cm ²
SPOT SIZE	15*26mm/15*40mm
REPETITION FREQUENCY	1-10Hz
TEMPERATURE OF SAPPHIRE	ICE PLUS, -42degree
INPUT POWER	110-120VAC, 50/60Hz
COOLING SYSTEM	Water+Air+Semiconductor

MACHINE PRINCIPLE




How does laser hair removal work?

Laser hair removal uses a laser that removes the hair follicle by targeting the melanin in the follicle. Melanin is the pigment that provides colour to the hair. The presence of melanin allows the laser to precisely target, heat and destroy the hair follicle.


This is clinic result of hair removal.

It takes 3-5 times to achieve permanent hair removal.

First treatment



Second treatment



Third treatment

