


CO₂ Fractional Laser Machine

Type Vertical

3 IN 1 Treatment System



Uniquely designed scanner head



7-joint light guide arm imported from South Korea



Laser power
70W



Focus spot diameter
20*20mm



Wavelength
10600nm





GLAMOURESSANCE

EMPOWERING BEAUTY. ENHANCING GLAMOUR.

ADVANTAGES



1. Laser power: **70W**



2. Laser source: RF (metal) tube laser imported from the United States; **5-10 years long life**





4. Cooling system: Air cooling



3. Light guide arm: **7-joint** light guide arm imported from South Korea




5. Wavelength: **10600nm**



6. Four models: **Fractional, High Energy Implantation, Surface Mode, Wrinkle Removal**



9. Scan graphics: Square, Circle, Diamond, Triangle, Ellipse, Straight line, Rectangle



7. Focus spot diameter: 20*20mm



8. Display: **11.6-inch** color touch screen



10. Scan mode: Random order, intermediate order, normal order

MACHINE FUNCTIONS



INFRARED AMING BEAM, FOR A PRECISE TREATMENT

FASTER SHOTS WITH LESS PAIN



Skin renewing and resurfacing



Skin tightening
wrinkle removal



Pore minimization



Sun damage recovery



Acne removal



Remove the pigment



Smooth scars and dermabrasion



Vaginal tightening,
Vulve Rejuvenation.




GLAMOURESSENCE
EMPOWERING BEAUTY. ENHANCING GLAMOUR.

LIGHT GUIDE ARM

South Korea imported 7 Articular Optical-arm.
Easily Operated And Greatly Reduces The Energy Loss



1. 360° rotation without obstacles



2. Uniquely designed scanner head More comfortable to hold



3. Imported from Korea, it has a longer service life and less energy loss.





GLAMOURESSENCE™

EMPOWER YOUR BEAUTY. UPLIFT YOUR GLAMOUR.

UPDATED PROBES


Variety of Scan Shapes And Handpieces Available



Korea imported high performance handpieces.

1 Fractional CO2 laser probe

For skin rejuvenation



Treatment graphic

Suitable for many different areas of treatment.



2 Normal CO2 laser probe


For surgery cutting



Safe Precision



No Downtime
Comfortable



A

B

C

3 Vaginal probes

For vaginal skin rejuvenation


A. Vulva Rejuvenation Probe

B. Treatment of Cervical Erosion Probe


C. Inner Vagina Rejuvenation Probe



1.Laser irradiation angle can be adjusted



2.Laser irradiation can be rotated 360 degrees



3.The area of the laser is compact edge contraction



4.Vaginal walls shrink and tightening


BEFORE AFTER




SCREEN DISPLAY

11.6-inch color touch screen


Mode Selection



Manual Mode



Smart Mode



MACHINE PARAMETERS

FRACTIONAL CO2 LASER SCAR REMOVER
 LASER RESURFACING



THE LIGHT GUIDE SYSTEM

7 Joint Imported Laser Arm

LASER SOURCE

RF Tube

FOCUS SPOT DIAMETER

20*20mm

PULSE WIDTH

0.1ms - 100ms

LASER POWER

70W

WAVELENGTH

10600nm

SCREEN

11.6 Inch